## Avansee™ Preload1P / Avansee™ Preload1P Clear

## Overview

Présentation / Überblick / Panoramica / Overzicht / Opis ogólny / Обзор / Vista general

1

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
**CHARACTERISTICS:**  
**Model YP2.2 (Yellow Type) / CP2.2 (Clear Type)**  
**OPTIC**  
 Material: Soft acrylic (UV-absorbing acrylic resin)  
 Colour: Yellow (YP2.2), Clear (CP2.2)  
 UV cutoff at 10% T: 409 nm (+20.0 dioptre lens of YP2.2), 397 nm (+20.0 dioptre lens of CP2.2)  
 Index of Refraction: 1.519 (35 °C)  
 Configuration: Biconvex  
 Power: +6.0 through +26.0 dioptre; +6.0 to +10.0 dioptre powers in +1.0 dioptre increments; +10.0 to +26.0 dioptre powers in +0.5 dioptre increments  
**HAPTICS**  
 Configuration: Modified-C  
 Material: Soft acrylic (UV-absorbing acrylic resin)  
 Colour: Yellow (YP2.2), Clear (CP2.2)

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
Figure 3 displays the typical spectral transmittance curves for YP2.2 / CP2.2 with a dioptric power of +20.0 together with the spectral transmittance curve\* for the phakic eye of a 53-year-old patient.

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




Insertion of the nozzle tip with the bevel down through an incision to just before the central pupillary area.

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







## PRESCRIPTION DEVICE

Caution: This device is restricted to sale by or on the order of an ophthalmologist.

## DEVICE DESCRIPTION

The Avansee™ Preload1P / Avansee™ Preload1P Clear from Kowa Company, Ltd. is a foldable posterior chamber aspheric intraocular lens (IOL) (Fig.2) that is preloaded in a single-use Kowa original injector (Fig.1). The optic of the lens and the modified-C haptics are made from a UV-absorbing hydrophobic soft acrylic material which, in the case of the Yellow type also contains proprietary blue-light filtering. The spectral transmittance of the IOL (Yellow type) closely replicates that of the natural crystalline lens (Fig.3). Prior to insertion, the optical portion of the lens is folded to allow for placement through a small incision. After surgical insertion into the eye, the lens gently unfolds to restore the optical performance.

## INTENDED USE

Kowa's Avansee™ Preload1P / Avansee™ Preload1P Clear IOL is intended to be positioned in the posterior chamber of the eye, replacing the natural crystalline lens. This allows the lens to function as a refractive medium in the visual correction of aphakia.

## INDICATIONS

Kowa's Avansee™ Preload1P / Avansee™ Preload1P Clear IOL is placed in a capsular bag and is designed for implantation after extracapsular cataract extraction or phacoemulsification of cataracts.

## WARNINGS

- Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting the Avansee™ Preload1P / Avansee™ Preload1P Clear lens in a child or a "special care patient" who has one or more of the following conditions. Before surgery Corneal endothelial damage / Glaucoma / Uveitis / Diabetic retinopathy / Retinal detachment / Congenital ocular anomalies / Choroidal haemorrhage / Shallow anterior chamber / Microphthalmos / Corneal dystrophy / Optic nerve atrophy / Ocular hypertension / Amydriasis / Amblyopia / History of keratoplasty / Iritis / Corneal anomalies / Macular degeneration / Retinal degeneration / Atopic disease / Pseudoexfoliation syndrome and zonular weakness / Zonular rupture and lens luxation (including lens subluxation) / Rubeosis iridis / Patients with intraoperative occurrences of any severe adverse event / or patients who, assessed by a surgeon, require special care for reasons such as an accompanying systemic or ocular disease.
- Prior to surgery, the risks and benefits associated with the implantation of the IOL should be clearly explained by the surgeon to the patient.
- Additional attention is required for "special care patients", including a post-operative follow up by an experienced ophthalmologist with adequate equipment, due to the higher risk of complications and/or insufficient vision recovery. Implantation of the IOL in children should be performed by an ophthalmologist who has sufficient knowledge and experience in paediatric care. This is particularly important when treating a child (under 2 years old) whose small eyeballs can make implantation and the handling of instruments difficult. The chances of IOL exchange are also higher due to the change of eye axial length with the child's growth. Therefore it is important to provide the parents of the young patient with informed consent information prior to any surgery and the IOL should be implanted to the young patients carefully after considering the benefit/risk ratio.
- For active uveitis and a child with uveitis, inflammation should be suppressed by medical treatment prior to intraocular lens implantation, as surgical invasion may cause aggravation of the uveitis or another complication.
- As with any surgical procedure, potential complications accompanying IOL implantation can occur. Adverse events and malfunctions may require discontinuation of the implantation, since they can lead to blindness, permanent problems with visual acuity and grave health hazard or may require IOL extraction and IOL exchange.
- Complications accompanying IOL implantation may include but are not limited to:
  - <Adverse events>
  - Corneal oedema / Keratitis (including corneal erosion) / Corneal endothelial damage / Acute corneal decompensation / Detachment of Descemet's membrane / Conjunctivitis and subconjunctival haemorrhage / Hyphaema / Hypopyon / Iris damage / Iritis (iridocyclitis) / Iris adhesion / Iris prolapse / Abnormal pupil (block, capture, deformation, dilatation, etc.) / Uveitis / Zonular rupture / Posterior capsular rupture / After cataract / Vitreous haemorrhage and opacity / Vitreous prolapse / Detachment, hole, tear, etc. of retinal tissue (including macula) / Retinal detachment / Choroidal detachment / Choroid haemorrhage / Macular oedema and degeneration / Expulsive haemorrhage / Endophthalmitis / Fibrin precipitation /

Secondary glaucoma / Intraocular pressure elevation (including transient elevated intraocular pressure and ocular hypertension) / Intraocular pressure lowering / Dysphotopsia / Deterioration of visual function (visual acuity and contrast sensitivity) / Error of predicted refractive power / Wound leak.

## &lt;Malfunctions&gt;

Optic damage (breakage, scratch, etc.) / Haptic damage (breakage, scratch, detachment, etc.) / Adhesion of foreign bodies on the lens surface / Lens surface reflection / Lens discoloration or pseudocoloration / Lens opacification (including glistening) / Lens luxation / Lens decentration / Lens dislocation / Lens jamming.

- Back-ups for the intraocular lens and the insertion device should be made available prior to surgery as a part of emergency preparations.
- DO NOT re-use the lens. For single use only.
- DO NOT re-sterilise. Re-use or re-sterilisation of any component may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury or illness. Re-use or re-sterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to the transmission of infectious disease(s) from one patient to another.
- DO NOT store the lens in direct sunlight or in any hot and humid places. Store the lens at room temperature (between 1 °C and 30 °C) and keep it dry.

## PRECAUTIONS

- Before using this device, please ensure that you have carefully read and understood this document in order to complete the procedure safely.
- This device should be placed in the posterior chamber. Implantation of the lens in the anterior chamber has not been evaluated for safety and efficacy.
- Before opening the Avansee™ Preload1P / Avansee™ Preload1P Clear packaging, confirm the description including the IOL model, lens power, and expiration date located on the label.
- The temperatures of the device and the ophthalmic viscosurgical device (OVD) during use should be kept between 21 °C and 25 °C.
- The inside of the aluminium packaging is not sterile. In a non-sterile environment, remove its contents including the primary (blister) package from the aluminium packaging.
- Should any damage or abnormality be found on the primary package after opening the IOL aluminium packaging, DO NOT use, as the device may no longer be sterile.
- After opening the primary package, the device should be handled aseptically.
- When removing the device from the primary package, ensure to be careful not to drop it.
- Before using the device, carefully examine all parts for damage, or any other abnormality. Should any abnormality be found, DO NOT use it.
- If the lens stage is dislodged before use, DO NOT use it.
- DO NOT remove the lens stage before the injection of the OVD.
- DO NOT open the cap on the injector body.
- If there are any abnormalities during the procedure, immediately stop using this device.
- When implantation has been stopped during the procedure due to an abnormality, DO NOT re-use the device and discard it.

## INSTRUCTIONS FOR USE

## &lt;Preparation for lens implantation&gt;

- Prior to implanting, examine the lens package for type, power and proper configuration.
- Open the primary (blister) package and remove the device in a sterile environment.
- Insert the injection needle for the OVD into the inlet located in the cap on the injector body, and fill the nozzle close to the dashed line shown in the diagram (Fig. 4) with the OVD in order to inject at least approx. 0.17 ml.

## &lt;Lens implantation&gt;

- With the Kowa mark portion of the injector body firmly in hand, hold the lens stage between your thumb and middle finger on both sides and support its fore part with your index finger (Fig. 5), and then remove the lens stage slowly away from the injector body towards the nozzle tip (Fig. 6).
- Push the plunger at a slow and constant rate; move the IOL forward and stop it at the point where the IOL optic is rolled and its edges make secure contact (Fig. 7). Once this has been done, immediately (within 20 seconds) place the IOL into the eye.
- Insert the nozzle tip with the bevel down through an incision to just before the central pupillary area (Fig. 8).
- While pushing the plunger ahead at a constant rate release the IOL inside the capsular bag.

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