

# Overview

Présentation / Überblick / Panoramica / Overzicht / Opis ogólny / Επισκόπηση / Vista general

# Specification

Spécifications / technische Beschreibung / Specifiche / Specificatie / Specyfikacja / Προδιαγραφές / Especificaciones

(a) Cylinder axis marks (at posterior surface)  
(a) Marque de l'axe cylindrique (sur la surface postérieure)  
(a) Markierungen der Zylinderachse (auf hinterer Oberfläche)  
(a) Segni asse del cilindro (nella superficie posteriore)  
(a) Merktelkens cilinders (op achteroppervlak)  
(a) Oznaczenie osi cylindra (na dolnej powierzchni)  
(a) Σημάδια του άξονα του κυλίνδρου (στην οπίσθια επιφάνεια)  
(a) Marcos del eje cilíndrico (en la superficie posterior)  
(b) Aspheric (anterior) surface  
(b) Surface asphérique (antérieure)  
(b) Asphärische (vordere) Oberfläche  
(b) Superficie (anterior) asferică  
(b) Asferisch (voor) oppervlak  
(b) Powierzchnia asferyczna (przednia)  
(b) Ασφαιρική (πρόσθια) επιφάνεια  
(b) Superficie asférica (anterior)  
(c) Toric (posterior) surface  
(c) Surface torique (postérieure)  
(c) Torische (hintere) Oberfläche  
(c) Superficie (posterior) torica  
(c) Torisch (achter) oppervlak  
(c) Powierzchnia toryczna (tylna)  
(c) Τορική (οπίσθια) επιφάνεια  
(c) Superficie tórica (posterior)  
Physical characteristics of the lens of YP-Tx / CP-Tx (x: 2, 3, 4, 5, 6, 7, 8 or 9)  
Caractéristiques physiques de la lentille de YP-Tx/CP-Tx (x: 2, 3, 4, 5, 6, 7, 8 or 9)  
Physikalische Eigenschaften der Linse von YP-Tx / CP-Tx (x: 2, 3, 4, 5, 6, 7, 8 or 9)  
Caratteristiche fisiche della lente di YP-Tx / CP-Tx (x: 2, 3, 4, 5, 6, 7, 8 or 9)  
Fysieke eigenschappen van de lens van de YP-Tx / CP-Tx (x: 2, 3, 4, 5, 6, 7, 8 or 9)  
Właściwości fizyczne soczewek YP-Tx / CP-Tx (x: 2, 3, 4, 5, 6, 7, 8 lub 9)  
Φυσικά χαρακτηριστικά του φακού του YP-Tx / CP-Tx (x: 2, 3, 4, 5, 6, 7, 8 ή 9)  
Características físicas de la lente de YP-Tx / CP-Tx (x: 2, 3, 4, 5, 6, 7, 8 o 9)

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







## PRESCRIPTION DEVICE

Caution: This device is restricted to sale by or on the order of an ophthalmologist.

## DEVICE DESCRIPTION

The Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear from Kowa Company, Ltd. is a foldable posterior chamber aspheric intraocular lens (IOL) (Fig.2) that is preloaded in a single-use Kowa original injector (Fig.1). The optic of the lens and the modified-C haptics are made from a UV-absorbing hydrophobic soft acrylic material which, in the case of the Yellow type also contains proprietary blue-light filtering. The spectral transmittance of the IOL (Yellow type) closely replicates that of the natural crystalline lens (Fig.3). Prior to insertion, the optical portion of the lens is folded to allow for placement through a small incision. After surgical insertion into the eye, the lens gently unfolds to restore the optical performance. The anterior surface of the lens is designed as aspheric and the cylinder power is added to posterior surface. Cylinder axis marks indicating flattest meridian are printed on the toric surface.

## INTENDED USE

Kowa's Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear IOL is intended to be positioned in the posterior chamber of the eye, replacing the natural crystalline lens. This allows the lens to function as a refractive medium in the visual correction of aphakia with corneal astigmatism.

## INDICATIONS

Kowa's Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear IOL is placed in a capsular bag and is designed for implantation after extracapsular cataract extraction or phacoemulsification of cataracts.

## WARNINGS

- Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting the Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear lens in a child or a "special care patient" who has one or more of the following conditions. Before surgery Corneal endothelial damage / Glaucoma / Uveitis / Diabetic retinopathy / Retinal detachment / Congenital ocular anomalies / Choroidal haemorrhage / Shallow anterior chamber / Microphthalmos / Corneal dystrophy / Optic nerve atrophy / Ocular hypertension / Amydriasis / Amblyopia / History of keratoplasty / Iritis / Corneal anomalies / Macular degeneration / Retinal degeneration / Atopic disease / Pseudoexfoliation syndrome and zonular weakness / Zonular rupture and lens luxation (including lens subluxation) / Rubeosis iridis / Patients with intraoperative occurrences of any severe adverse event / or patients who, assessed by a surgeon, require special care for reasons such as an accompanying systemic or ophthalmic disease.
- Prior to surgery, the risks and benefits associated with the implantation of the IOL should be clearly explained by the surgeon to the patient.
- Additional attention is required for "special care patients", including a post-operative follow up by an experienced ophthalmologist with adequate equipment, due to the higher risk of complications and/or insufficient vision recovery. Implantation of the IOL in children should be performed by an ophthalmologist who has sufficient knowledge and experience in paediatric care. This is particularly important when treating a child (under 2 years old) whose small eyeballs can make implantation and the handling of instruments difficult. The chances of IOL exchange are also higher due to the change of eye axial length with the child's growth. Therefore it is important to provide the parents of the young patient with informed consent information prior in any surgery and the IOL should be implanted to the young patients carefully after considering the benefit/risk ratio.
- For active uveitis and a child with uveitis, inflammation should be suppressed by medical treatment prior to intraocular lens implantation, as surgical invasion may cause aggravation of the uveitis or another complication.
- Rotation of the lens away from its intended axis may reduce the astigmatic correction. Misalignment of 30 degrees or greater may deteriorate postoperative astigmatism. If reoperation for axis correction is required, it should be performed as early as possible before the lens adheres to the capsule bag.
- As with any surgical procedure, potential complications accompanying IOL implantation can occur. Adverse events and malfunctions may require discontinuation of the implantation, since they can lead to blindness, permanent problems with visual acuity and grave health hazard or may require IOL extraction and IOL exchange.
- Complications accompanying IOL implantation may include but are not limited to:
  - <Adverse events>
  - Corneal oedema / Keratitis (including corneal erosion) / Corneal endothelial damage / Acute corneal decompensation / Detachment of Descemet's membrane / Conjunctivitis and subconjunctival haemorrhage / Hypohaema / Hypopyon / Iris damage / Iritis (iridocyclitis) / Iris adhesion / Iris prolapse / Abnormal pupil (block,

capture, deformation, dilatation, etc.) / Uveitis / Zonular rupture / Posterior capsular rupture / After cataract / Vitreous haemorrhage and opacity / Vitreous prolapse / Detachment, hole, tear, etc. of retinal tissue (including macula) / Retinal detachment / Choroidal detachment / Choroid haemorrhage / Macular oedema and degeneration / Expulsive haemorrhage / Endophthalmitis / Fibrin precipitation / Secondary glaucoma / Intraocular pressure elevation (including transient elevated intraocular pressure and ocular hypertension) / Intraocular pressure lowering / Dysphotopsia / Deterioration of visual function (visual acuity and contrast sensitivity) / Error of predicted refractive power / Wound leak.

## <Malfunctions>

Optic damage (breakage, scratch, etc.) / Haptic damage (breakage, scratch, detachment, etc.) / Adhesion of foreign bodies on the lens surface / Lens surface reflection / Lens discoloration or pseudocolouration / Lens opacification (including glistening) / Lens luxation / Lens decentration / Lens dislocation / Lens jamming.

- Back-ups for the intraocular lens and the insertion device should be made available prior to surgery as a part of emergency preparations.
- DO NOT re-use the lens. For single use only.
- DO NOT re-sterilise. Re-use or re-sterilisation of any component may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury or illness. Re-use or re-sterilisation may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to the transmission of infectious disease(s) from one patient to another.
- DO NOT store the lens in direct sunlight or in any hot and humid places. Store the lens at 30 °C or lower and keep it dry.

## PRECAUTIONS

- Before using this device, please ensure that you have carefully read and understood this document in order to complete the procedure safely.
- DO NOT use this lens for the patients with capsular rupture or with zonular rupture, or with large posterior capsulotomy.
- Appropriate IOL model and placement of axis should be determined. The location and size of the incision may affect the amount and axis of residual corneal astigmatism. Kowa provides a calculation tool ([www.avansee.eu/](http://www.avansee.eu/)) for surgeons and recommends its use.
- This device should be placed in the posterior chamber. Implantation of the lens in the anterior chamber has not been evaluated for safety and efficacy.
- The lens should not be placed in the ciliary sulcus.
- Before opening the Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear packaging, confirm the description including the IOL model, dioptric and cylinder lens power, and expiration date located on the label.
- The temperatures of the device and the ophthalmic viscosurgical device (OVD) during use should be kept between 21 °C and 25 °C.
- The inside of the aluminium packaging is not sterile. In a non-sterile environment, remove its contents including the primary (blister) package from the aluminium packaging.
- Should any damage or abnormality be found on the primary package after opening the aluminium packaging, DO NOT use, as the device may no longer be sterile.
- After opening the primary package, the device should be handled aseptically.
- When removing the device from the primary package, ensure to be careful not to drop it.
- Before using the device, carefully examine all parts for damage, or any other abnormality. Should any abnormality be found, DO NOT use it.
- If the lens stage is dislodged before use, DO NOT use it.
- DO NOT remove the lens stage before the injection of the OVD.
- DO NOT open the cap on the injector body.
- If there are any abnormalities during the procedure, immediately stop using this device.
- Except for aborting lens insertion, do not interrupt insertion after the lens starts coming out from the nozzle, and keep pushing the plunger until trailing haptic is completely extruded from the nozzle.
- When implantation has been stopped during the procedure due to an abnormality, DO NOT re-use the device and discard it.

## INSTRUCTIONS FOR USE

<Preparation for lens implantation>

- Prior to implanting, examine the lens package for type, power and proper configuration.
- Open the primary (blister) package and remove the device in a sterile environment.
- Insert the injection needle for the OVD into the inlet located in the cap on the injector body, and fill the nozzle close to the dashed line shown in the diagram (Fig. 4) with the OVD in order to inject at least approx. 0.17 ml.

<Lens implantation>

- With the *Kowa* mark portion of the injector body firmly in hand, hold the lens stage between your thumb and middle finger on both sides and support its fore part with your index finger (Fig. 5), and then remove the lens stage slowly away from the injector body straight (Fig. 6).
  - Push the plunger at a slow and constant rate; move the IOL forward and stop it at the point where the IOL optic is rolled and its edges make secure contact (Fig. 7). Once this has been done, immediately (within 20 seconds) place the IOL into the eye.
  - Insert the nozzle tip with the inlet (*Kowa* mark) upward through an incision to just before the central pupillary area (Fig. 8).
  - While pushing the plunger ahead at a constant rate release the IOL inside the capsular bag.
  - Push until the trailing haptic is completely extruded and the IOL is released, and pull the injector body out of the anterior chamber. Then, stabilise the IOL inside the capsular bag using a hook or other tools.
  - Care should be taken to remove the OVD from the eye and align the cylinder axis marks with the intended axis at the end of surgery.
- Cautions in handling:**
- A sufficient amount of OVD should be injected to fill the nozzle and cover the entire lens optic by inserting the injection needle into the inlet located in the cap on the injector body (to the level of at least 0.17 ml from the inlet). Insufficient OVD may damage the IOL because of inadequate lubrication inside the nozzle.
  - DO NOT leave the injector with the OVD filled inside the nozzle. Immediately start the next procedure.
  - DO NOT advance the plunger before removing the lens stage.
  - When removing the lens stage, DO NOT twist it laterally (Fig. 9).
  - Upon removing the lens stage, confirm that there is no abnormality in the IOL. Should a haptic of the IOL become deformed and protrude out of the hole in the back of the injector body as shown in Fig. 10, DO NOT use it.
  - After removing the lens stage, DO NOT fill it with additional OVD.
  - Push the plunger slowly when it makes contact with the lens.
  - DO NOT push the plunger rapidly or intermittently. Push the plunger forward at a smooth constant rate.
  - DO NOT pull back the plunger once the plunger has been pushed forward.
  - Keep the inlet (*Kowa* mark) upward while pushing the plunger ahead at a constant rate to release the IOL inside the capsular bag.
  - If you feel strong resistance or any abnormal feeling during advancement of the plunger, stop using this device.
  - During advancement of the plunger, stop using this device if a trailing haptic is out of the groove along the plunger, or if :
    - the leading haptic should twist and extend forward (Fig. 11);
    - the leading haptic should bend and stretch out (Fig. 12);
    - the trailing haptic should extend out of the nozzle (Fig. 13);
    - the plunger should pass above or under the lens optic or bend the optic irregularly (Fig. 14); or
    - the plunger should go too far toward the right or left side (Fig. 15).
  - The IOL should be moved forward until the lens optic is rolled and its edges make secure contact. If you stop advancement short of this

location, the plunger may get into and be caught in the lens optic and the lens may be damaged when pushing and inserting the IOL afterward.

- After the IOL is moved into the near tip of the nozzle, move on to the next procedure promptly (within 20 seconds) and place the lens into the eye. Failure of prompt implantation may cause the lens to become stuck inside the nozzle, possibly resulting in damage to the IOL or improper implantation.
- During the IOL implantation procedure, DO NOT push the lens to the extent of being exposed from the nozzle tip (Fig. 16).
- The IOL should be released slowly after the nozzle tip has been located just before the centre of the pupillary area.
- All residual OVD on and under the lens should be carefully removed, otherwise it may cause lens rotation and misalignment.

**CALCULATION OF LENS POWER**

The surgeon should determine preoperatively the spherical equivalent and cylindrical power of the lens to be implanted. The A-constant and the cylinder powers on IOL plane below are presented as a reference value for lens power calculations. When calculating the exact lens power, it is recommended that calculations should be performed individually based on the equipment used and the operating surgeon's experience.

Model	A-constant (A-mode)
YP-Tx / CP-Tx	118.6

**Cylinder powers (D):**

Model : YP-/CP-	T2	T3	T4	T5	T6	T7	T8	T9
IOL plane (labeled)	0.75	1.50	2.25	3.00	3.75	4.50	5.25	6.00
Corneal plane	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0

Lens power calculation methods are described in the following references.

Sanders, D.R., Retzlaff, J., and Kraff, M.C., "Comparison of SRK II formula and other second generation formulae", *Journal of Cataract and Refractive Surgery*, Vol. 14, pp.136-141, 1988.

Retzlaff, J.A., Sanders, D.R., and Kraff, M.C., "Development of the SRK/T intraocular lens implant power calculation formula", *Journal of Cataract and Refractive Surgery*, Vol. 16, pp.333-340, 1990; ERRATA, Vol. 16, pp.528, 1990.

**INTRAOCULAR LENS PATIENT CARD**

The Intraocular Lens Patient Card included in the outer packaging should be completed and given to the patient. The patient must provide this card when seeking any future medical diagnosis.

**REPORT OF INCIDENT**

In the event of any incident that may reasonably be regarded as device related, please contact your local distributor or Authorized Representative (See below).

**HOW SUPPLIED**

The Avansee™ Preload1P Toric / Avansee™ Preload1P Toric Clear is supplied sterile in a primary (blister) package for single use only. The primary (blister) package is sterilized with ethylene oxide gas and the package should be opened only under sterile conditions. DO NOT re-sterilise.

**HOW DISCARDED**

After the lens implantation, please discard the device as medical waste in accordance with local laws and regulations.

**EXPIRATION DATE (sterility expiration date)**

The sterility expiration date is clearly indicated on the outer package. The date is expressed as four digits for the year and two digits for the month.

**RETURN/EXCHANGE POLICY**

Regarding lens return or exchange, please contact with your local distributor or Authorized Representative (See below).

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Symbol Explanation:

Symbols	Description	Symbols	Description	Symbols	Description
	Do not re-use / For single use only		Caution		Manufacturer
	Do not re-sterilise		Keep away from sunlight		Authorized representative in the European Community
	Do not use if package is damaged		Keep dry		D
	Sterilised using ethylene oxide		Temperature limitation		Overall diameter
	Non-sterile		Serial number		Optic diameter (Body diameter)
	Consult instructions for use		Use by (YYYY-MM : year-month)		Information website for patients
	Patient Name or patient ID		Date of implantation		Name and Address of the implanting healthcare institution
	Medical device		UDI as AIDC format		Importer

**Manufacturer:**  
*Kowa Company, Ltd.*  
4-14, 3-Chome, Nihonbashi-Honcho, Chuo-ku, Tokyo 103-8433, JAPAN

**Authorized Representative in the European Community:**  
*Kowa Pharmaceutical Europe GmbH*  
Görreshof 151, 53347 Alfthor, Germany  
Phone: 00-800-2826-7333  
Fax: +49 2222 95891 99  
[www.kowapharmaceuticals.eu](http://www.kowapharmaceuticals.eu)  
[www.avansee.eu](http://www.avansee.eu)

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